

OSG Newsletter October 2016

In this Newsletter:

- OSG releases a new version of **RT Software Suite**! Check the video for an overview of the key features.
- In November OSG participates in many **events**: the 1st Dutch Sleep Conference, the KNF Technician's Day in Groningen and the NVLA Conference in Ede. See you soon!

Curious about the New Features of our Latest Version?

Check out the video!

RT Software Suite 3.2



New features



Release 3.2 also brings the following improvements:

- We have optimized the import of **Embla & Embletta** signals for frequent use;
- Parameters for Respiratory analysis and Arousal detection have been finetuned;
- You can filter Requests more easily on their **Location**;
- High Frequency Oscillations (**HFO**) **filters** are available in the Derivation Editor;
- We've optimized RT Software Suite for **ThinApp** virtualization;
- We support two new Video Formats for Ambulatory Video: MTS and MPE.

Log in to the OSG website <http://www.osg.be> for the Application Notes and the complete list of improvements!

Interested in a Software Upgrade? Contact OSG!