

OSG

List of PSG supplies for Morpheus








Table of Contents

Morpheus accessories.....	3
Thermocan cannula/kabel.....	6
Nonin sensors.....	7
SleepSense inductive bands and sensors	9
Medcare locks	11
Technomed electrodes / cables	12
EEG cup electrodes	12
Cables	13
Electrode linkers	14
Disposable EEG cup electrodes	15
Disposable EMG electrodes.....	16
Cannulas (disposables)	17




Morpheus accessories



Product code	Description
0602/ECG EL TPCLB	ECG Cable with connectors Touch Proof - Clip, length <u>150 cm</u> , Black 
0602/ECG EL TPCLR	ECG Cable with connectors Touch Proof - Clip, length <u>150 cm</u> , Red
0539/1041200-Morph	ECG cables: Reusable snap leads, DIN 42 802 connector, <u>76 cm</u> , per set of 2 pieces (black and white) 
0602/BSP OXY	Oxymeter to be used with MORPHEUS Recorder, with cable  
0602/BQ USB	USB Interface Card for MORPHEUS Serie Headbox 
0602/HBC INT LTM R	Cable to Connect Interfaces to MORPHEUS or to Quick Disconnect TERM LTM, New Model

List of PSG supplies for Morpheus




	
0602/TERM LTM	Terminal to connect with MORPHEUS, with Quick Disconnect 
0602/BAG + BELT MORPH	Bag and belt for Morpheus recorders  
0602/BQ STAND MORPH	MORPHEUS Headbox Stand, composed of BASE, STAND and Headbox Support SUPP MORPH 
0539/P1223	Pro-Tech airflow sensor 1 channel pediatric 2 thermocouples nasal + 1 mouth 
0539/P1214	Pro-Tech airflow sensor, 1 channel pediatric small, nasal bridgemount

		
0539/P1273	Pro-Tech Airflow sensor 1 element	
0464/0514	Braebon cTherm Airflow Sensor (Cannula Compatible Airflow Sensor)	


Thermocan cannula/kabel

0454/15807	Thermocan 7-Ft Adult Thermal Cannula (<u>per stuk</u>) 
0454/1475	Adult Thermocan Interface Cable(Thermocouple) / Safety DIN 
0454/15813	Thermocan 7-Ft Pediatric Thermal Cannula (<u>per stuk</u>) 
0454/1475P	Pediatric Thermocan Interface Cable(Thermocouple) / Safety DIN 

Nonin sensors

Product code	Description
0433/OX.8000J	<p>Nonin Adult Flex Sensor Includes 8000J sensor, attachment tape (25) and 1 meter cable</p> 
0433/OX.8000J-3	<p>Nonin Adult Flex Sensor Includes 8000J sensor, attachment tape (25) and 3 meter cable</p>
0433/OX.8008J	<p>Nonin Infant Flex Sensor Includes 8008J sensor, Hydrogel tape strips (25) and 1 meter cable</p> 
0433/OX.8001J	<p>Nonin Neonatal Flex Sensor Includes 8001J sensor, Hydrogel tape strips (25) and 1 meter cable</p> 
0433/OX.8000JFW	Nonin flexi wrap adult attachment tape for Nonin 8000J sensor (25pcs)
0433/OX.8008JFW	Nonin flexi wrap pediatric attachment tape for Nonin 8008J sensor (25 pcs)
0433/OX.8000Q	Nonin ear clip sensor

List of PSG supplies for Morpheus

		
--	---	--



SleepSense inductive bands and sensors

Product code	Description
0454/9002-L40	Multiple use inductive band - small (2/pack) 
0454/9002-L60	Multiple use inductive band - medium (2/pack) 
0454/9002-L90	Multiple use inductive band - large (2/pack) 
0454/9002-L120	Multiple use inductive band - x-large (2/pack) 
0454/9010	Reusable Inductive Effort Belt (Blue) (<u>per piece</u>) (1 year warranty) 
0454/9010-KIT	Reusable Inductive Effort Belt Kit (2 Belts + 9101 + 9102)
0454/9141	RIP Interface Cable 1.5 ft (45cm) - Abdomen / Key Connector 
0454/9142	RIP Interface Cable 1.5 ft (45cm) - Chest / Key Connector

List of PSG supplies for Morpheus


	
0454/9151	RIP Interface Cable 8 ft (2,4m) - Abdomen / Key Connector
0454/9152	RIP Interface Cable 8 ft (2.4m) - Chest / Key Connector
0454/9101	Inductive Interface Cable - Abdomen / Safety Connectors 
0454/9102	Inductive Interface Cable - Chest / Safety Connectors
0454/1250	Piezo Snore Sensor / Safety DIN connectors 
0454/14807	Luer to O2 adapter (pack of 5) 
0454/14808	Male to male Luer adapter(pack of 5) 

Medcare locks

Product code	Description
0611/1421020_2	Medcare locks (XactTrace Lock – Abdomen + Thorax) 
0611/1421020	XactTrace Lock – Abdomen
0611/1421022	XactTrace Lock – Thorax
0611/1421000	XactTrace belt, 20m 

Technomed electrodes / cables

EEG cup electrodes



EEG cup electrodes have been designed to withstand long-term use by medical technicians. The precision stamped electrodes are available in two diameters: the pediatric, 6mm (0.24") or the adult, 10mm (0.39") size. They are stamped in two materials, 99.9% pure silver or tin. The silver electrodes are then covered with a heavy gold-plate, silver chlorided or kept silver. The flexible silicon coated wires are standard in multi colors.


<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C12-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C12-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C12-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector
0234/C14-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C14-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C22-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C22-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C22-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

List of PSG supplies for Morpheus

0234/C24-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C24-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C24-834	per set of 12 pieces	(6 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

<u>CODE</u>	<u>PACKING</u>	<u>DESCRIPTION</u>
0234/C32-434	per set of 12 pieces	(10 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C32-634	per set of 12 pieces	(10 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C32-834	per set of 12 pieces	(10 mm cup); 2 m silicon cable multi color with DIN 42 802 connector
0234/C34-434	per set of 12 pieces	(6 mm cup); 1 m silicon cable multi color with DIN 42 802 connector
0234/C34-634	per set of 12 pieces	(6 mm cup); 1,5 m silicon cable multi color with DIN 42 802 connector
0234/C34-834	per set of 12 pieces	(6 mm cup); 2 m silicon cable multi color with DIN 42 802 connector

Cables



Product code	Description
0539/1041200	ECG cables: Reusable snap leads 5 colors, DIN 42 802 connector, 76 cm, per set of 5 pieces 
0539/1075233	ECG cables: Reusable snap leads 5 colors, DIN 42 802 connector, 183 cm, per set of 5 pieces

		
0539/1075400	EMG cables: Reusable snap leads, DIN 42 802 connector, 183 cm, per set of 10 pieces	
0539/1075401	EMG cables: Reusable snap leads, DIN 42 802 connector, 243 cm, per set of 10 pieces	
0234/D93-832	Reusable Extension cable 2m multicolor male-female DIN42802 connector 6 pc / package	

Electrode linkers

Product code	Description
0234/ET30-101	Electrode linker red
0234/ET30-201	Electrode linker black

List of PSG supplies for Morpheus

	
0234/ET21-008RD	Double female adapter red, in: 2 mm pin connector, out: DIN 42 802 connector (per 25 pcs)
	
0234/ET22-008	Female/male adapters, 2mm pin connector (per 25 pcs)


Disposable EEG cup electrodes



Product code	Packaging	Description
0234/C52-438	5x5 pieces	DISPOSABLE CUP AG/CL ELEKTRODES 10MM DIN42802 – 1meter
0234/C52-838	5x5 pieces	DISPOSABLE CUP AG/CL ELEKTRODES 10MM DIN42802 – 2meter

Disposable EMG electrodes



Product code	Packaging	Description
0684/72010-K/10	Per 420	Neuroline 720 EMG elektrode, 22x30mm met kabel 100cm en K-connector (1,5mm)
0684/72020-K/C/12	Per 480	Neuroline 720 EMG elektrode, 22x30mm met kabel 200cm en K-connector (1,5mm)
0234/TE/K50430-002	Per 40	DISPOSABLE ADHESIVE SURFACE ELEKTRODES 

Cannulas (disposables)



Product code	Description
0539/P1259	ProTech adult sensor cannulas (60 per case)
0539/P1295	Pro-Tech adult nasal/oral cannulas (30 per case)
0539/P1263	Pro-Tech pediatric nasal sensor cannulas (60 per case)
0539/P1267	Pro-Tech pediatric small nasal sensor cannulas (60 per case)
0539/P1324	Cannula adapter kit (3 per set) <div data-bbox="459 1189 742 1357" data-label="Image"> </div>
0539/P1107	Cannula Titration Filter (50 pack) <div data-bbox="488 1473 651 1675" data-label="Image"> </div>